



# SURFACE™



## Technical Specification Sheet

SurFACE™ FABRIC SWAB PRE-MOISTENED WITH NEUTRALIZING BUFFER



|                    |  |
|--------------------|--|
| Item No.           | 10006188   |
| Description        | Fabric Swab (Swab Cloth according to ISO 18593) in Twirl-Tie™ bag, individually packed & dosed with 10 ml of NB  |
| Intended Use       | For environmental surfaces to collect and ship microbial samples   |
| Components         | Fabric Swab (Swab Cloth) 50% wood pulp/50% polyester, size: 28x35cm, 4-fold, weight: 75gsm, pre-moistened with 10 ml buffer, 55 oz. (1.62 L) Twirl-Tie™ Bag, 7.5" x 12" (19 x 30 cm)   |
| Buffer             | NEUTRALIZING BUFFER<br><br>Lecithin 7 g/L,<br>Polysorbate 80 (Tween 80) 5 g/L<br>Potassium phosphate di-basic 8.5 g/L<br>Sodium Bisulfite 1.25 g/L<br>Sodium Phosphate di-basic 3 g/L<br>Sodium Thioglycolate 1g/L<br>Sodium Thiosulfate 6 g/L |
| Sterilization      | Sterilized by irradiation at dosage of 25 KGy  |
| Shelf Life         | 12 months from date of production  |
| Storage Conditions | 4-25 °C; protect from light, heat and moisture   |
| Pack Size          | 200 in a box   |
| Quality Control    | Recovery/inhibitory of microorganisms<br>Sterility   |
| Conformity         | Product is conform with ISO 18593:2018 "Microbiology of the food chain - Horizontal methods for surface sampling"  |

Date of Creation: 05. May 2020